

Optimizing the Performance of a LaBr₃:Ce, Sr Based Neutron and Gamma Dual Detector

Kan Yang*, Peter R. Menge

Saint-Gobain Crystals, HIRAM, OH, USA *kan.yang@saint-gobain.com

Introduction

LaBr₃:Ce scintillators, when co-doped with a small amount of strontium and other aliovalent elements, are shown to have significantly improved light output and energy resolution [1, 2]. Our previous study shows strontium and calcium co-doping also enhances the crystal's ability to discriminate between heavy charged particles and gamma radiations via pulse shape discrimination (PSD). Complete suppression of internal alpha radiation background can be achieved in Ca or Sr co-doped LaBr₃:Ce [3]. Moreover, the enhanced alpha-gamma PSD makes it possible to use the co-doped LaBr₃:Ce crystal for dual neutron-gamma detection when coupled with a neutron converter such as ⁶LiF, which is also an excellent reflector.

Pulse Shape Discrimination




Fig. 1 Averaged PMT pulses of (a) Ce only and (b) Sr co-doped LaBr₃:Ce when excited by a particles or γ-rays

- ❖ Sr co-doped LaBr₃:Ce shows more prominent pulse shape differences between gamma and alpha pulses than standard Ce only LaBr₃.

Detector Configuration




Fig. 2 Schematic drawing of the prototype detector




Fig. 3 Photo of the assembled detector (w/o PMT)

$n + {}^6\text{Li} \rightarrow t (2.75 \text{ MeV}) + \alpha (2.05 \text{ MeV})$

- ❖ 93.5% enriched ⁶LiF powder is used as thermal neutron converter.
- ❖ Tritons (2.75 MeV) from ⁶Li(n,t)α reaction are used as the signal for neutron detection.
- ❖ MCNPX simulations suggest the optimal ⁶LiF layer thickness is 60 μm.
- ❖ Detectors with different ⁶LiF thickness (30 μm, 60 μm, 100 μm, and 200 μm) were tested.




Fig. 4 Averaged neutron energy deposition per incident neutron (simulated)

Thermal Neutron Detection

Experimental Setup A ²⁵²Cf source (12.8 ng) was placed 5" from the detector. The source was shielded with 2.25" thick Pb and 2" thick high density polyethylene (HDPE) to attenuate the gammas and moderate the neutrons, respectively. Thermal neutron flux through the detector was determined by MCNPX.




Fig. 5 Gamma-neutron PSD scatter intensity plot of LaBr₃:Ce, Sr with (a) 30 μm (b) 60 μm (c) 100 μm (d) 200 μm thick ⁶LiF conversion layer

Table I Neutron detection properties of ⁶LiF surrounded LaBr₃:Ce,Sr detectors (LLD = 1.5 MeV)

| ⁶ LiF Thickness | Thermal Neutron Detection ($E_n < 0.2 \text{ eV}$) | | | PSD FoM | Gamma Rejection Ratio* |
|----------------------------|--|----------------------|--------------------------------|---------|------------------------|
| | Capture Rate | Detection Efficiency | Det. Eff. per unit area | | |
| 30 μm | 8.1% | 2.9% | 0.60% /layer*inch ² | 1.19 | 5.2×10^{-5} |
| 60 μm | 15.2% | 3.1% | 0.63% /layer*inch ² | 1.17 | 7.8×10^{-5} |
| 100 μm | 23.3% | 2.8% | 0.56% /layer*inch ² | 1.23 | 2.2×10^{-6} |
| 200 μm | 39.8% | 2.7% | 0.55% /layer*inch ² | 1.23 | 1.9×10^{-6} |

* Assuming 95% of neutrons are counted

- ❖ As expected, thermal neutron detection efficiency is highest when the thickness of ⁶LiF is 60 μm.
- ❖ A fraction of captured neutrons generate detectable signals depending on the thickness of ⁶LiF.
- ❖ Further increasing the thickness of ⁶LiF will not increase neutron detection efficiency.




Fig. 6 Pulse height spectra for tritons from ⁶Li(n,t)α reaction, showing increasing straggling with ⁶LiF thickness




Fig. 7 Pulse shape spectra of 60 μm ⁶LiF coated detector at different lower energy thresholds

- ❖ Reducing ⁶LiF thickness reduces energy straggling of tritons in ⁶LiF.
- ❖ Minimal change in triton energy spectrum if ⁶LiF thickness > 100 μm.
- ❖ Good PSD when lower energy threshold > 1.5 MeV.




Fig. 8 a) Thermal neutron detection efficiencies ($E_n < 0.2 \text{ eV}$), b) PSD Figure of Merit, and c) gamma rejection ratio at different lower energy thresholds

Summary and Outlook

- ❖ A dual gamma-neutron detector can be created by surrounding a Sr co-doped LaBr₃:Ce crystal with ⁶LiF powder.
- ❖ A fast PMT is not required to realize good PSD in Sr co-doped LaBr₃:Ce. Main pulse shape differences are in the slow decay component.
- ❖ Multiple ⁶LiF layers or large surface-area-to-volume ratio are desired for efficient neutron detection.
- ❖ Replacing the usual reflector with ⁶LiF powder can add neutron detection capability to LaBr₃:Ce,Sr detectors.

REFERENCES

- [1] K. Yang, P.R. Menge, J. Buzniak, V. Ouspenski, 2012 IEEE-NSS-MIC Conference record, N1-135, pp. 308-311 (2012).
- [2] M. S. Alekhin, J. T. M. de Haas, I. V. Khodyuk, K. W. Krämer, P. R. Menge, V. Ouspenski and P. Dorenbos, Appl. Phys. Lett. 102, 161915 (2013).
- [3] K. Yang, P.R.Menge, V. Ouspenski, 2014 IEEE-NSS-MIC Conference record, N39-2 (2014)