

Improved light collimation for scintillators crystals using a photonic crystal


Why a photonic crystal (PhC) ?

We aim to collimate light extracted from a scintillator crystal, in order to improve coupling when standoff between the scintillator and the detector is required (eg. high energy radiography with an electron accelerator, IRM-PET, reactor core imaging and human PET imaging).

Currently light extraction from scintillators reaches more than 80% thanks to:

- wrapping of high reflective Teflon film
- roughened surfaces in order to reflect light with nearly random angles

It's hard to further increase the quantity of the extracted light ... however there is substantial room for improving angular distribution.



To change Fresnel law at the interface we propose a Photonic crystal, which is a 2D grating with period close to the extracted wavelength


PhC effect - Numerical calculations

We consider a cylindrical scintillator crystal of diameter 50 mm, height 50 mm and optical index $n=1,85$. Using the ray tracing software, LighTools, we calculated the extracted power from one of the circular faces of the scintillator collected by a far-field detector


PhC were simulated with a Rigorous Couple Wave Analysis (RCWA) method

PhC parameters are:

- Period (a)
- Pillar width (d)
- Pillar height (h)



We calculated the fraction of extracted optical power in the solid angle as a function of vertex semi-angle Θ with respect to a flat surface (gain). We studied the gain stability as a function of the PhC main parameters:




A PhC with index 1.85, $a=400$, $d=200$ nm, $h=200$ nm can increase small angle light extraction up to 50%.


PhC fabrication

The proposed fabrication technique to realize high index photonic crystals is based on embossing a SiO_2 matrix doped with TiO_2 nanoparticles. Embossing is a relatively fast and low-cost process compared to electron beam lithography and this technique paves the way to large scale fabrication.

By changing the concentration of the TiO_2 nanoparticles we can finely tune the optical index of the layer (see first panel Fig.3). We selected a concentration of 60% to obtain an index of 1.85




We demonstrated the experimental fabrication of a high index photonic crystal by embossing a solgel containing TiO_2 nanoparticles




We measure a shrinkage of about 30%. This has to be taken into account during the mold design

Picture of the PhC diffracting green light (The scintillator is covered with a black mask)



Preliminary measurement

Light extracted from the PhC in a solid angle of 1° is five times larger than a smooth surface



Conclusions and perspectives

We designed and fabricated a high index PhC on a LYSO scintillator crystal

- We will soon optically characterize our sample in order to verify simulations reliability by comparing them with experimental data and therefore experimentally demonstrate the collimation effect.
- In the future we aim at emboss a larger surface of many square centimeters

[1] M. Kronberger et al. "Improving Light Extraction from Heavy Inorganic Scintillators by Photonic Crystals", IEEE Transactions on Nuclear Science (Vol 57, Issue 5)
 [2] A. Knapitsch "Photonic crystals: A novel approach to enhance the light output of scintillation based detectors" Nuclear Instruments and Methods in Phy. Res. A
 [3] Pawel Modrzynski et al. "Light extraction from scintillating crystals enhanced by photonic crystal structures patterned by focused ion beam" IEEE transaction on nuclear science
 [4] Patent WO2015197947 "Luminescent material with a textured photonic layer"